


CleanMark Autoclave Indicator Labels

High Performance Labels

CleanMark at a glance

Many demanding applications found in the Life Sciences and Semiconductor industries call for performance requirements beyond the limits of standard off-the-shelf labels. At CleanMark, we can provide the most appropriate and compatible labels for various individual product and operational process needs.

Key Features



ISO 5

Cleanroom processed



121 °C

Activation temperature



Type 1

11140-1 compliant
process indicator



Superior

Moisture durability



CLEANMARK
HIGH PERFORMANCE LABELS

503-968-8303

www.cleanmarklabels.com


6700 SW Bradbury CT, Portland, OR 97224

THE CHALLENGE

Steam sterilization via high-pressure autoclaving is critical for eradicating biological pathogens in healthcare and medical device settings. Operators must ensure their items reach the proper temperature for a specified time to ensure sterilization. During the process, labels must withstand the extreme conditions of the chamber.

THE SOLUTION

CleanMark autoclave labels are constructed to remain intact and legible while undergoing steam sterilization. Our autoclave compatible labels play a critical role in assisting operators identify what is being sterilized, when it was processed, and how it should be safely handled. Additionally, CleanMark autoclave labels may feature a color changing indicator to aid operators verify whether the conditions to sufficiently inactivate biological pathogens have been met. CleanMark autoclave labels are customizable for a wide range of specific requirements and working conditions.



**Highly Uniform
Color Contrast**



**Direct to Page
Indicator Printing**



**Available in
white to blue**



**Water-based
non-toxic ink**

CUSTOMIZATION

Our autoclave labels are fully customizable to satisfy a wide range of technical requirements. We will walk you through all of the autoclave label considerations to develop you a complete label solution.